

รายละเอียดและคุณสมบัติอุปกรณ์

- | | |
|--|---------------|
| 1. Mini-Mill PULVERISETTE23 | Quantity x2 |
| <ul style="list-style-type: none">● Final finess: 5 μm● Oscillations per minute: 900-3000 at 9 mm amplitude● Power 220V 50 Hz | |
| 2. Zirconium oxide grinding bowl | Quantity x2 |
| <ul style="list-style-type: none">● volume: 10 ml | |
| 3. Zirconium oxide grinding bowl | Quantity x2 |
| <ul style="list-style-type: none">● volume: 15 ml | |
| 4. Zirconium oxide grinding ball | Quantity x100 |
| <ul style="list-style-type: none">● diameter: 5 mm | |
| 5. Zirconium oxide grinding ball | Quantity x10 |
| <ul style="list-style-type: none">● diameter: 10 mm | |

Remark: Item 2-5 can be use with item 1.